

Continenence



MED ITALIA SERVICE

info@meditaliaservice.com
meditaliaservicesrls@pec.it

+39 3314605453
 Via delle barozze 6 A
 00040 - Rocca di Papa (Rm)

www.meditaliaservice.com
shop.meditaliaservice.com

Partita iva: 16193371008
 Numero REA: RM - 1640791
 SDI: M5UXCR1

Dual channel NMES device for Incontinence

Basic Pelvic Toner to promote Continenence.
 3 Custom programmes plus 9 Preset programmes.

Pelvic Floor programmes:

- Pelvic Floor pain relief
- Urge, Stress, Frequency Incontinence
- Lack of Sensitivity
- Pelvic Floor Work Out
- Endurance

Lock mode function to measure the home compliance (mA and time used).

Features

STIM/TENS Channels

9

Pre-Set Programs

3

Custom Programs

Battery Level

PP3

Power Supply

Home Compliance

Waveform

Belt Clip

Accessories supplied with the device

- Dual lead wire x 2
- 9V PP3 battery
- Soft bag
- Operating Manual
- Gift Box
- Optional:
- Vaginal Probe

Technical Data

Output current	0-90 mA (with 500Ω load)	Output modes	Continuos, Work/Rest
Waveform	Asymmetrical, rectangular	Power Supply	PP3 9V battery
Pulse Frequency	2- 100 Hz		
Phase Duration	50-450 µs		
Weight	90 g		