


Anal Ultrasound Probe



Visual Value...
..Innovative 360° Transrectal Ultrasound



Fecal Incontinence Imaging

Features and Benefits:

- High resolution USB-enabled 360° 12MHz endocavity probe - the smallest in the world
- Anal canal visualization capability for diagnostics and surgical planning
- Integrated fecal incontinence diagnostics
- Ultrasound Video Capture Capability
- Self-contained hardware for simple cleaning and service
- One-Click image marking with customizable annotations
- Easy to use angle and length measurements
- Intuitive hot keys for image quality adjustment
- Extract images for use in MS PowerPoint, Word and PDF
- Optional HL7 connectivity



Code: 99-5935

SPECIFICATIONS:


Frequency:	12MHz
Scan Depth:	3cm Max
Focal Point:	1.6cm
Probe Length:	31.5cm
Probe Diameter:	2.1cm
Weight:	0.4kg
Connectivity:	USB 2.0
Depth Marker:	2-12cm

SOFTWARE:

Measurements:	Dist. (mm), Angle (°s)
Image Capture:	Frames & Cine Loop
Cine Loop Length:	128, 64 & 28 Frames



Case Study #1: Normal



Case Study #2: Sphincter Defect



MED ITALIA SERVICE



info@meditaliaservice.com
meditaliaservicesrls@pec.it



www.meditaliaservice.com
shop.meditaliaservice.com



+39 3314605453
Via delle barozze 6 A
00040 - Rocca di Papa (Rm)

Partita iva: 16193371008
Numero REA: RM - 1640791
SDI: M5UXCR1



INTERSON