



Rehabilitation
Pain Management

Verity Medical Ltd.
*Designer and Manufacturer
of NeuroTrac devices*

All our products are manufactured in compliance with the Medical Device Directive
MDD 93/42/EEC

Developed and Manufactured in the UK



MyoPlus2 / MyoPlus4



The EMG and ETS tool

Two channel EMG tool for Biofeedback training and neuro-muscular assessment with easy to read backlit LCD display. Advanced ETS used for neuro re-education: the bar graph shows the muscular effort and the stimulation is triggered automatically to support the effort.

Pre-set and customised programmes

Two or four channel STIM devices with 3 built-in preset clinical modes: Continenence, Sports, Rehab. Each of these modes has a set of programmes with up to 5 time phases. Each phase is EMG or STIM or ETS.

The device also includes 4 custom clinical modes: EMG Custom, STIM Custom, ETS Custom and P+A Custom.

Wireless (Bluetooth) connection

Up to 10 metres of wireless (Bluetooth) connection for NeuroTrac Software.

Features



Accessories supplied with the device

- Dual lead wire x 2 or 4 (depending on model)
- EMG reference wire
- Self adhesive electrodes (50 x 50 mm)
- 1.5 V AA batteries x4
- Carrying Case
- Stand
- Rubber sleeve with belt clip
- Operating Manual
- Gift Box

Technical Data

EMG Range	0.2 - 2000 μ V RMS
Sensitivity	0.1 μ V RMS
Output current	0-90 mA (with 500 Ω load)
Waveform	Symmetrical, rectangular
Pulse Frequency	2- 100 Hz
Phase Duration	50-450 μ s

Weight	160 g
Power Supply	4 x 1.5 V AA batteries



NeuroTrac[®]



MED ITALIA SERVICE



info@meditaliaservice.com
meditaliaservicesrls@pec.it



www.meditaliaservice.com
shop.meditaliaservice.com



+39 3314605453
Via delle barozze 6 A
00040 - Rocca di Papa (Rm)

Partita iva: 16193371008
Numero REA: RM - 1640791
SDI: M5UXCR1



Verity Medical Ltd
Unit 7, Upper Slackstead Farm
Farley Lane, Braishfield
Hampshire SO51 0QL
United Kingdom