



BIOLITO

MYOLITO

CONTINENCE

PELVIFINE

PELVITONE

EVOSTIM-UG

n. channels	2	2	2	2	2	2
stim currents	cont work / rest (WR) modulated - burst	cont work / rest (WR) modulated - burst	cont work / rest (WR)	cont work / rest (WR)	cont work / rest (WR) modulated	cont work / rest (WR)
tens currents	han	han	sì	sì	sì	sì
free programs	1	1	3+1 (cont)	3	3	10
n ° sessions free pr.	1/5	1/5	1	1/5	1/5	1/10
accessories	4 electrodes 40x40 4 electrodes 40x80	4 electrodes 40x40 4 electrodes 40x80	no	no	no	MINIMA vaginal probe
number of programs	5	5 + 7	9	11	11	5x10
Urge incontinence	1	1	1	1	1	1
Stress/effort incontinence Faecal incontinence	2	2	2	2	2	1
Mixed incontinence	1	1	no	no	no	1
Pain reliever	1	1	1	1	1	1
Frequency / Urgency	no	no	2	3	3	no
Lack of sensitivity	no	no	1	1	1	no
Pelvic floor (toning)	no	no	1	1	1	1
Relaxation	no	no	no	1	1	no
Endurance-Resistance	no	no	1	1	1	no
Nerve stimulation	no	4	no	no	no	no
Muscle stimulation	no	3	no	no	no	no