

## KegelSmart™

EAN: 7350022 278004


Get stronger, healthier pelvic floor muscles with the easiest Kegel coach ever. KegelSmart™ uses touch sensor technology to register your pelvic floor strength and then sets the perfect exercise routine for you. All you have to do is: contract when it vibrates and rest when it stops. Developed with gynecologists and made with the smoothest medical grade silicone, KegelSmart gives you an expert pelvic floor workout in less than 5 minutes a day!



Features	<ul style="list-style-type: none"> <li>Registers pelvic strength</li> <li>Automatically sets exercise level</li> <li>Adjusts the routine as you improve</li> <li>Tracks your progress through 5 levels</li> </ul>
Benefits	<p>KegelSmart™ helps women of all ages to:</p> <ul style="list-style-type: none"> <li>Strengthen their pelvic floor</li> <li>Prepare for a healthy pregnancy</li> <li>Recover vaginal tightness after childbirth</li> <li>Enhance sensations during intimacy</li> </ul>
Includes	<ul style="list-style-type: none"> <li>1 KegelSmart™</li> <li>Storage pouch</li> <li>User manual</li> <li>Requires: 1 AAA battery</li> </ul>

### Easiest Kegel Routine Ever

- 1 Turn it on and insert
- 2 **Contract** when it vibrates  
**Rest** when it stops
- 3 The device turns off automatically once the routine is over



### Product Specifications

Product with Packaging

Dimensions: 115 x 48 x 140 mm / 4.53 x 1.89 x 5.51 in.

Weight: 147 g / 5.19 oz.

### Shipment Specifications

	Dimensions	Weight	Product Count
<b>Shipping Box</b>	585 x 355 x 260 mm 23.03 x 13.98 x 10.24 in.	8.63 kg / 19.02 lbs	48 units

#### Intimina Headquarters EU Office

LELO AB, Karlavägen 41,  
114 31 Stockholm, Sweden

#### US Office

LELO Inc., 5799 Fontanoso Way  
San Jose, CA 95138, USA

#### Asia Office

LELO, Room 1701-1710, 75 South  
Wan Ping Road Shanghai 200032, China